

Appendix A: Cross-Network PBMC Processing Worksheet v7.0

Note: The fields in this worksheet must be filled out by hand, using a pen.

Specimen Processing Laboratory:		Protocol:		
Participant ID (PTID/PID):		Visit Number:	Visit Type:	
Collection Date:		Collection Time:		
Processing Start Date:		Processing Start Time:	Processed By (Initials):	
Reagents	Manufacturer	Lot Number		Expiration Date
DMSO				
FBS				
WDR: HBSS or PBS (circle one)				
Cell Separation Tube (frit)				
Density Gradient Media				
	Volume in mL (record as X.Y)			Expiration Date
CPS	CPS	DMSO	FBS	1 working day (<18hrs)
Data to be Captured During Processing				Sample
Sample tube type (circle one or record "other" tube type)				ACD / HEP / EDT Other: _____
Blood condition (circle one or more; add comments on reverse as needed)				SAT/ HEM / CLT
Measured usable whole blood volume (to the nearest 0.1mL)				mL
Indicate processing method (circle one)				CSTFB / overlay / underlay
Counting Method: Name of specific instrument or manual count (record in field to right)				
Counting re-suspension volume of HBSS (or other WDR) (V) (record as X.Y)				mL
Cell count average concentration (C)				x 10 ⁶ cells/mL
Total cell number (T) = C x V				x 10 ⁶ cells
Calculate cell yield/mL of whole blood. (QC check)= (T/Usable Whole Blood Volume)				x 10 ⁶ cells/mL
If T/A * ≥ N1 ; CPS re-suspension vol (V1) = A If T/A < N1 ; Calculate estimated CPS re-suspension vol. (V1)=(T/N1x10⁶ cells/mL)(1mL)				mL
Calculate final CPS re-suspension volume (V_f), (V1 rounded DOWN to nearest whole (X.0) mL)				mL
Calculate actual number of cells per vial. N2 = (T/V_f) x V2 ; (V2 =1 mL) Note: Do not store more than 50 million cells per vial				x 10 ⁶ cells/vial
Print and QC LDMS Label content/barcodes (initials of person (s) performing QC)				
Frozen Date and Time (ddMMMyyyy /HH:MM) (Explain in comments section if not within 4 hours of processing start time)				
Number of Cryovials actually frozen Note: Should be equal to final CPS re-suspension volume for 1mL aliquots (V_f) and ≤ (A)				
Complete remaining LDMS entries including total cell count & frozen time (Initials)				

*Note: **A** = The maximum number of aliquots required according to the protocol-specific Laboratory Processing Chart (LPC).
Do not store more than this number of aliquots.

Appendix A: Cross-Network PBMC Processing Worksheet v7.0 (Page 2 of 2)

Note: The fields in this worksheet must be filled out by hand, using a pen.

Specimen Processing Laboratory:

PTID/PID:

Transfer of Cryovials to Freezer Storage Box	
Person who transferred cryovials to storage box locations assigned by LDMS	
Date (ddMMMyyyy)/time cryovials were transferred from controlled-rate freezing device to storage box. (Sample must be maintained at -80°C during transfer)	
Initial (Primary) Review (Initials/Date)	
Final (Secondary) Review (Initials/Date)	

Hemocytometer Counts	Total Count	Viable Cells	Non-Viable
Square #1 (cells/mm ²)			
Square #2 (cells/mm ²)			
Square #3 (cells/mm ²)			
Square #4 (cells/mm ²)			
Average Cell Count per Square (cells/mm ²)			
PBMC Dilution Factor (1:DF**)			
Hemocytometer Factor for cells/mL	10 ⁴	10 ⁴	10 ⁴
Cell count concentration (C) = (Average Cells/mm ²)(DF)(10 ⁴); convert to 10 ⁶ cells/mL	Not applicable	x 10 ⁶ cells/mL	Not applicable
% viability = (Viable cells 4 squares/total cells 4 squares) (100)	Not applicable		Not applicable

**Note: Dilution Factor (DF) = (parts cells + parts dilution fluid)/ parts cells

Automated Cell Counts (10 ³ /μl=10 ⁶ /mL)	Count #1
Cell Count (C) as cells x 10 ⁶ /mL	
PBMC Dilution Factor (1:DF***)	
Cell Concentration = (C)(DF)	x 10 ⁶ cells/mL
% viability (if applicable)	

***Note: Dilutions for automated counters are extremely rare. If performing direct counts, enter a 1 in the DF box and complete the column.

Comments, protocol deviations, and additional information not captured elsewhere in this worksheet:

Note: (A) = The maximum number of aliquots required according to the protocol-specific Laboratory Processing Chart (LPC). Do not store more than this number of aliquots.

Cross-Network PBMC Processing Standard Operating Procedure

Appendix A: Cross-Network PBMC Processing Worksheet v7.0 Ex: N1 = 20 mil cells; A = 5 aliquots

Note: The fields in this worksheet must be filled out by hand, using a pen.

Specimen Processing Laboratory: **Lab 398**

Protocol: **313**

Participant ID (PTID/PID): **123-456789**

Visit Number: **2.0**

Visit Type: **vst**

Collection Date: **08AUG2024**

Collection Time: **08:00**

Processing Start Date: **08AUG2024**

Processing Start Time:
08:45

Processed By (Initials): **CN**

Reagents	Manufacturer	Lot Number	Expiration Date
DMSO	Sigma	RNBM0548	18JAN2025
FBS	Peak	13G1212	18AUG2025
WDR: HBSS or PBS (circle one)	Gibco	2660057	30APR2026
Cell Separation Tube (frit)	Greiner	E220337Q	14MAR2027
Density Gradient Media	Cytiva	1Q345061	31AUG2026
Volume in mL (record as X.Y)			Expiration Date
CPS Prepared 19AUG2024 08:30 CN	CPS	DMSO	FBS
	9.0	0.9	8.1
			1 working day (<18hrs)
Data to be Captured During Processing			Sample
Sample tube type (circle one or record "other" tube type)			<input checked="" type="radio"/> CD / HEP / EDT
Blood condition (circle one or more; add comments on reverse as needed)			Other: _____
			<input checked="" type="radio"/> AT / HEM / CLT
Measured usable whole blood volume (to the nearest 0.1mL)			86.3 mL
Indicate processing method (circle one)			<input checked="" type="radio"/> CSTFB / <input type="radio"/> overlay / underlay
Counting Method: Name of specific instrument or manual count (record in field to right)			Manual Count
Counting re-suspension volume of HBSS (or other WDR) (V) (record as X.Y)			17.0 mL
Cell count average concentration (C)			7.2 x 10 ⁶ cells/mL
Total cell number (T) = C x V			122.4 x 10 ⁶ cells
Calculate cell yield/mL of whole blood. (QC check)= (T/Usable Whole Blood Volume)			1.4 x 10 ⁶ cells/mL
If T/A ≥ N1; CPS re-suspension vol (V1) = A			5.0 mL
If T/A < N1; Calculate estimated CPS re-suspension vol. (V1)=(T/N1x10 ⁶ cells/mL)(1mL)			5.0 mL
Calculate final CPS re-suspension volume (V _f), (V ₁ rounded DOWN to nearest whole (X.0) mL)			5.0 mL
Calculate actual number of cells per vial. N2 = (T/V _f) x V2; (V2=1 mL)			24.4 x 10 ⁶ cells/vial
Note: Do not store more than 50 million cells per vial			
Print and QC LDMS Label content/barcodes (initials of person (s) performing QC)			CN
Frozen Date and Time (ddMMMyyyy/HH:MM) (Explain in comments section if not within 4 hours of processing start time)			08AUG2024
Number of Cryovials actually frozen			5
Note: Should be equal to final CPS re-suspension volume for 1mL aliquots (V _f) and ≤ (A)			
Complete remaining LDMS entries including total cell count & frozen time (Initials)			CN

***Note:** A = The maximum number of aliquots required according to the protocol-specific Laboratory Processing Chart (LPC). Do not store more than this number of aliquots.

Example #1:

N1 = 20x10⁶ cells/mL

A = 5 aliquots

Calculations:

CPS re-suspension volume (V1)

$$122.4/5 = 24.4 > 20$$

Thus, T/A ≥ N1

$$(V1) = A$$

Actual number of cells per vial (N2)

$$122.4/5 \times 1 =$$

$$24.4 \times 10^6 \text{ cells/vial}$$

Appendix A: Cross-Network PBMC Processing Worksheet v7.0 Ex: N1 = 10 mil cells; A = 5 aliquots

Note: The fields in this worksheet must be filled out by hand, using a pen.

Specimen Processing Laboratory: **Lab 398**

Protocol: **313**

Participant ID (PTID/PID): 123-456789		Visit Number: 2.0		Visit Type: vst
Collection Date: 08AUG2024		Collection Time: 08:00		
Processing Start Date: 08AUG2024		Processing Start Time: 08:45		Processed By (Initials): CN
Reagents	Manufacturer	Lot Number	Expiration Date	
DMSO	Sigma	RNBM0548	18JAN2025	
FBS	Peak	13G1212	18AUG2025	
WDR: HBSS or PBS (circle one)	Gibco	2660057	30APR2026	
Cell Separation Tube (frit)	Greiner	E220337Q	14MAR2027	
Density Gradient Media	Cytiva	1Q345061	31AUG2026	
	Volume in mL (record as X.Y)			Expiration Date
CPS Prepared 19AUG2024 08:30 CN	CPS	DMSO	FBS	1 working day (<18hrs)
	9.0	0.9	8.1	
Data to be Captured During Processing				Sample
Sample tube type (circle one or record "other" tube type)				ACD / HEP / EDT Other: _____
Blood condition (circle one or more; add comments on reverse as needed)				SAT / HEM / CLT
Measured usable whole blood volume (to the nearest 0.1mL)				46.3 mL
Indicate processing method (circle one)				CSTFB / overlay / underlay
Counting Method: Name of specific instrument or manual count (record in field to right)				Manual Count
Counting re-suspension volume of HBSS (or other WDR) (V) (record as X.Y)				9.0 mL
Cell count average concentration (C)				4.2 x 10⁶ cells/mL
Total cell number (T) = C x V				37.8 x 10⁶ cells
Calculate cell yield/mL of whole blood. (QC check)= (T/Usable Whole Blood Volume)				0.8 x 10⁶ cells/mL
If T/A * ≥ N1; CPS re-suspension vol (V1) = A				3.7 mL
If T/A < N1; Calculate estimated CPS re-suspension vol. (V1)=(T/N1x10 ⁶ cells/mL)(1mL)				3.0 mL
Calculate final CPS re-suspension volume (V _f), (V1 rounded DOWN to nearest whole (X.0) mL)				3.0 mL
Calculate actual number of cells per vial. N2 = (T/V _f) x V2; (V2=1 mL)				12.6 x 10⁶ cells/vial
Note: Do not store more than 50 million cells per vial				
Print and QC LDMS Label content/barcodes (initials of person (s) performing QC)				CN
Frozen Date and Time (ddMMMyyyy /HH:MM) (Explain in comments section if not within 4 hours of processing start time)				08AUG2024
Number of Cryovials actually frozen				3
Note: Should be equal to final CPS re-suspension volume for 1mL aliquots (V _f) and ≤ (A)				
Complete remaining LDMS entries including total cell count & frozen time (Initials)				CN

**Note:* A = The maximum number of aliquots required according to the protocol-specific Laboratory Processing Chart (LPC). Do not store more than this number of aliquots.

Example #2:

N1 = 10x10⁶

cells/mL

A = 5 aliquots

Calculations:

CPS re-suspension volume (V1)

$$37.8/5 = 7.56 < 10$$

Thus, T/A < N1

$$(V1) = 37.8/10 \times 10^6 \text{ cells/mL} (1\text{mL})$$

Actual number of cells per vial (N2)

$$37.8/3 = 12.6 \times 10^6 \text{ cells/vial}$$